

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RPTN Polyclonal Antibody

多克隆抗体

产品货号: PA05732

第三版

描述:RPTN (Repetin) is a Protein Coding gene. Diseases associated with RPTN include Autosomal Recessive Congenital Ichthyosis. Among its related pathways are Keratinization and Developmental Biology. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is FLG2.

抗原:Synthesized peptide derived from the C-terminal region of human Repetin

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expression is scattered in the normal epidermis but strong in the acrosyringium, the inner hair root sheath and in the filiform papilli of the tongue.