

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RFX3 Polyclonal Antibody

多克隆抗体

产品货号: PA05739

第三版

描述: This gene is a member of the regulatory factor X gene family, which encodes transcription factors that contain a highly-conserved winged helix DNA binding domain. The protein encoded by this gene is structurally related to regulatory factors X1, X2, X4, and X5. It is a transcriptional activator that can bind DNA as a monomer or as a heterodimer with other RFX family members. Multiple transcript variants encoding different isoforms have been described for this gene. RFX3 (Regulatory Factor X3) is a Protein Coding gene. Diseases associated with RFX3 include Testicular Lymphoma and Pineal Gland Cancer. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and transcription regulatory region DNA binding. An important paralog of this gene is RFX2.

抗原: Synthesized peptide derived from the C-terminal region of human RFX3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: