

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RHOD Polyclonal Antibody

多克隆抗体

产品货号: PA05746

第三版

**描述:** Ras homolog, or Rho, proteins interact with protein kinases and may serve as targets for activated GTPase. They play a critical role in muscle differentiation. The protein encoded by this gene binds GTP and is a member of the small GTPase superfamily. It is involved in endosome dynamics and reorganization of the actin cytoskeleton, and it may coordinate membrane transport with the function of the cytoskeleton. Two transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthesized peptide derived from the Internal region of human Rho D

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Heart, placenta, liver, skeletal muscle, and pancreas and, with weaker intensity, in several other tissues.