

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RHOBTB3 Polyclonal Antibody

多克隆抗体

产品货号: PA05748

第三版

**描述:** RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351). RHOBTB3 (Rho Related BTB Domain Containing 3) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and Retrograde transport at the Trans-Golgi-Network. GO annotations related to this gene include GTP binding and Rab GTPase binding.

**抗原:** Synthesized peptide derived from the Internal region of human RHOBTB3

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.