



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APBB1IP Polyclonal Antibody

多克隆抗体

产品货号: PA05752

第三版

**描述:**APBB1IP (Amyloid Beta Precursor Protein Binding Family B Member 1 Interacting Protein) is a Protein Coding gene. Among its related pathways are RET signaling and Innate Immune System. An important paralog of this gene is RAPH1.

**抗原:**Synthesized peptide derived from the Internal region of human RIAM

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Widely expressed with high expression in thymus, spleen, lymph node, bone marrow and peripheral leukocytes.