

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey RPL10L Polyclonal Antibody

多克隆抗体

产品货号: PA05753

第三版

描述: This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. This gene is intronless. RPL10L (Ribosomal Protein L10 Like) is a Protein Coding gene. Among its related pathways are Viral mRNA Translation and rRNA processing in the nucleus and cytosol. GO annotations related to this gene include structural constituent of ribosome. An important paralog of this gene is RPL10.

抗原: Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L10L

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Testis specific.