

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RIN3 Polyclonal Antibody

多克隆抗体

产品货号: PA05779

第三版

描述:Summary: This protein encoded by this gene is a member of the RIN family of Ras interaction-interference proteins, which are binding partners to the RAB5 small GTPases. The protein functions as a guanine nucleotide exchange for RAB5B and RAB31. Alternative splicing results in multiple transcript variants. RIN3 (Ras And Rab Interactor 3) is a Protein Coding gene. Diseases associated with RIN3 include Paget's Disease Of Bone. Among its related pathways are Vesicle-mediated transport and RAB GEFs exchange GTP for GDP on RABs. GO annotations related to this gene include GTPase activator activity and Rab guanyl-nucleotide exchange factor activity. An important paralog of this gene is RIN2.

抗原:Synthesized peptide derived from the Internal region of human RIN3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Widely expressed.