
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MPRIP Polyclonal Antibody

多克隆抗体

产品货号: PA05780

第三版

描述:MPRIP (Myosin Phosphatase Rho Interacting Protein) is a Protein Coding gene. Among its related pathways are HIV Life Cycle and Oncogenic MAPK signaling. GO annotations related to this gene include actin binding. An important paralog of this gene is TRIOBP.

抗原:Synthesized peptide derived from the Internal region of human M-RIP.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: