

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZC3H7B Polyclonal Antibody

多克隆抗体

产品货号: PA05789

第三版

描述: This gene encodes a protein that contains a tetratricopeptide repeat domain. The encoded protein also interacts with the rotavirus non-structural protein NSP3. ZC3H7B (Zinc Finger CCCH-Type Containing 7B) is a Protein Coding gene. GO annotations related to this gene include poly(A) RNA binding. An important paralog of this gene is ZC3H7A.

抗原: Synthesized peptide derived from the C-terminal region of human RoXaN

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: