

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse RASL10A Polyclonal Antibody

多克隆抗体

产品货号: PA05792

第三版

**描述:**RASL10A (RAS Like Family 10 Member A) is a Protein Coding gene. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is RASL10B.

**抗原:**Synthesized peptide derived from the Internal region of human RRP22

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expression appears to be strictly limited to the central nervous system.