

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RPS6KA6 Polyclonal Antibody

多克隆抗体

产品货号: PA05797

第三版

描述: This gene encodes a member of ribosomal S6 kinase family, serine-threonine protein kinases which are regulated by growth factors. The encoded protein may be distinct from other members of this family, however, as studies suggest it is not growth factor dependent and may not participate in the same signaling pathways.

抗原: Synthesized peptide derived from the C-terminal region of human Rsk-4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: