

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human S100Z Polyclonal Antibody

多克隆抗体

产品货号: PA05801

第三版

描述:Members of the S100 protein family contain 2 calcium-binding EF-hands and exhibit cell-type specific expression patterns. For additional background information on S100 proteins, see MIM 114085.S100Z (S100 Calcium Binding Protein Z) is a Protein Coding gene. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is S100A1.

抗原:Synthesized peptide derived from the N-terminal region of human S-100Z

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: