

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DLGAP1 Polyclonal Antibody

多克隆抗体

产品货号: PA05807

第三版

描述:DLGAP1 (DLG Associated Protein 1) is a Protein Coding gene. Diseases associated with DLGAP1 include Retinitis Pigmentosa 55 and Obsessive-Compulsive Disorder. Among its related pathways are Protein-protein interactions at synapses and Neuroscience. An important paralog of this gene is DLGAP2.

抗原:Synthesized peptide derived from the C-terminal region of human SAPAP1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in brain.