

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SCAF1 Polyclonal Antibody

多克隆抗体

产品货号: PA05809

第三版

**描述:**SCAF1 (SR-Related CTD Associated Factor 1) is a Protein Coding gene. GO annotations related to this gene include RNA binding and protein domain specific binding. An important paralog of this gene is PHRF1.

**抗原:**Synthesized peptide derived from the N-terminal region of human SCAF1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Ubiquitous. Highly expressed in fetal brain and liver, poorly expressed in salivary gland, heart, skin and ovary. Expressed in colorectal carcinomas and ovarian cancers. Overexpressed in colorectal carcinomas as compared to normal colonic mucosa.