

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINB4 Polyclonal Antibody

多克隆抗体

产品货号: PA05811

第三版

描述:The protein encoded by this gene is a member of the serpin family of serine protease inhibitors. The encoded protein is highly expressed in many tumor cells and can inactivate granzyme M, an enzyme that kills tumor cells. This protein, along with serpin B3, can be processed into smaller fragments that aggregate to form an autoantigen in psoriasis, probably by causing chronic inflammation. SERPINB4 (Serp Family B Member 4) is a Protein Coding gene. Diseases associated with SERPINB4 include Squamous Cell Carcinoma and Chromosome 18Q Deletion Syndrome. Among its related pathways are Amoebiasis. GO annotations related to this gene include enzyme binding and protease binding. An important paralog of this gene is SERPINB3.

抗原:Synthesized peptide derived from the Internal region of human SCCA2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Squamous cells.