

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human JTB Polyclonal Antibody

多克隆抗体

产品货号: PA05813

第三版

描述: Required for normal cytokinesis during mitosis. Plays a role in the regulation of cell proliferation. May be a component of the chromosomal passenger complex (CPC), a complex that acts as a key regulator of mitosis. The CPC complex has essential functions at the centromere in ensuring correct chromosome alignment and segregation and is required for chromatin-induced microtubule stabilization and spindle assembly. Increases AURKB activity. Inhibits apoptosis induced by TGFB1 (By similarity). Overexpression induces swelling of mitochondria and reduces mitochondrial membrane potential.

抗原: Synthesized peptide derived from the N-terminal region of human JTB

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous. Expressed in all normal human tissues studied but overexpressed or underexpressed in many of their malignant counterparts.