

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SCTR Polyclonal Antibody

多克隆抗体

产品货号: PA05817

第三版

描述:The protein encoded by this gene is a G protein-coupled receptor and belongs to the glucagon-VIP-secretin receptor family. It binds secretin which is the most potent regulator of pancreatic bicarbonate, electrolyte and volume secretion. Secretin and its receptor are suggested to be involved in pancreatic cancer and autism.

抗原:Synthesized peptide derived from the Internal region of human Secretin Receptor

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: