

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SEPT3 Polyclonal Antibody

多克隆抗体

产品货号: PA05828

第三版

描述: This gene belongs to the septin family of GTPases. Members of this family are required for cytokinesis. Expression is upregulated by retinoic acid in a human teratocarcinoma cell line. The specific function of this gene has not been determined. Alternative splicing of this gene results in two transcript variants encoding different isoforms. SEPT3 (Septin 3) is a Protein Coding gene. Diseases associated with SEPT3 include Teratocarcinoma and Sexual Disorder. Among its related pathways are Bacterial invasion of epithelial cells. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is SEPT9.

抗原: Synthesized peptide derived from the Internal region of human Septin 3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Brain-specific.