

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SEPT6 Polyclonal Antibody

多克隆抗体

产品货号: PA05829

第三版

**描述:** This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. SEPT6 (Septin 6) is a Protein Coding gene. Diseases associated with SEPT6 include Sexual Disorder. Among its related pathways are Bacterial invasion of epithelial cells. GO annotations related to this gene include GTP binding. An important paralog of this gene is SEPT11.

**抗原:** Synthesized peptide derived from the Internal region of human Septin 6

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Widely expressed.