



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINB12 Polyclonal Antibody

多克隆抗体

产品货号: PA05834

第三版

描述:SERPINB12 (Serp Family B Member 12) is a Protein Coding gene. Among its related pathways are Innate Immune System. GO annotations related to this gene include enzyme binding and serine-type endopeptidase inhibitor activity. An important paralog of this gene is SERPINB4.

抗原:Synthesized peptide derived from the Internal region of human Serpin B12

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in many tissues, including brain, bone marrow, lymph node, heart, lung, liver, pancreas, testis, ovary, and intestine.