

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SFRS14 Polyclonal Antibody

多克隆抗体

产品货号: PA05835

第三版

描述:This gene encodes a member of the arginine/serine-rich family of splicing factors. The encoded protein functions in mRNA processing. Alternatively spliced transcript variants have been described.SUGP2 (SURP And G-Patch Domain Containing 2) is a Protein Coding gene. Among its related pathways are mRNA Splicing - Major Pathway. GO annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is SUGP1.

抗原:Synthesized peptide derived from the Internal region of human SFRS14

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in adult testis, and in fetal brain and kidney.