

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SFRS15 Polyclonal Antibody

多克隆抗体

产品货号: PA05836

第三版

描述: This gene likely encodes a member of the arginine/serine-rich splicing factor family. A similar protein in Rat appears to bind the large subunit of RNA polymerase II and provide a link between transcription and pre-mRNA splicing. Alternatively spliced transcript variants have been described. SCAF4 (SR-Related CTD Associated Factor 4) is a Protein Coding gene. Diseases associated with SCAF4 include Hay-Wells Syndrome and Wells Syndrome. GO annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is SCAF8.

抗原: Synthesized peptide derived from the Internal region of human SFRS15

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: