

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SIAH1/2 Polyclonal Antibody

多克隆抗体

产品货号: PA05850

第三版

描述:SIAH1 (Siah E3 Ubiquitin Protein Ligase 1) is a Protein Coding gene. Diseases associated with SIAH1 include Parkinson Disease, Late-Onset. Among its related pathways are Proteolysis Role of Parkin in the Ubiquitin-Proteasomal Pathway and Innate Immune System. GO annotations related to this gene include identical protein binding and ubiquitin-protein transferase activity. An important paralog of this gene is SIAH2. SIAH2 (Siah E3 Ubiquitin Protein Ligase 2) is a Protein Coding gene. Among its related pathways are Proteolysis Role of Parkin in the Ubiquitin-Proteasomal Pathway and Innate Immune System. GO annotations related to this gene include ligase activity and transcription corepressor activity. An important paralog of this gene is SIAH1.

抗原:Synthesized peptide derived from the Internal region of human Siah-1/2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:SIAH1 is Widely expressed at a low level. Down-regulated in advanced hepatocellular carcinomas.SIAH2 isWidely expressed at low level.