

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SIX5 Polyclonal Antibody

多克隆抗体

产品货号: PA05858

第三版

描述:The protein encoded by this gene is a homeodomain-containing transcription factor that appears to function in the regulation of organogenesis. This gene is located downstream of the dystrophia myotonica-protein kinase gene. Mutations in this gene are a cause of branchiootorenal syndrome type 2.SIX5 (SIX Homeobox 5) is a Protein Coding gene. Diseases associated with SIX5 include Branchiootorenal Syndrome 2 and Branchiootorenal Syndrome. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is SIX4.

抗原:Synthesized peptide derived from the Internal region of human Six5

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in adult but not in fetal eyes. Found in corneal epithelium and endothelium, lens epithelium, ciliary body epithelia, cellular layers of the retina and the sclera.