

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC17A2 Polyclonal Antibody

多克隆抗体

产品货号: PA05860

第三版

描述:SLC17A2 (Solute Carrier Family 17 Member 2) is a Protein Coding gene. Among its related pathways are Osteoblast Signaling. GO annotations related to this gene include sodium:phosphate symporter activity. An important paralog of this gene is SLC17A4.

抗原:Synthesized peptide derived from the Internal region of human SLC17A2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in the small intestine, kidney, spleen and testis. Not detected in fetal brain, bone marrow, and mammary gland.