

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SLK Polyclonal Antibody

多克隆抗体

产品货号: PA05863

第三版

描述:SLK (STE20 Like Kinase) is a Protein Coding gene. Diseases associated with SLK include Amyotrophic Lateral Sclerosis 9. Among its related pathways are Apoptosis and Autophagy and Adhesion. GO annotations related to this gene include protein homodimerization activity and transferase activity, transferring phosphorus-containing groups. An important paralog of this gene is STK10.

抗原:Synthesized peptide derived from the C-terminal region of human SLK

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed. Highest expression is found in heart and in skeletal muscle.