

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PLK2 Polyclonal Antibody

多克隆抗体

产品货号: PA05883

第三版

描述:The protein encoded by this gene is a member of the polo family of serine/threonine protein kinases that have a role in normal cell division. This gene is most abundantly expressed in testis, spleen and fetal tissues, and its expression is inducible by serum, suggesting that it may also play an important role in cells undergoing rapid cell division. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. PLK2 (Polo Like Kinase 2) is a Protein Coding gene. Diseases associated with PLK2 include Kummell's Disease. Among its related pathways are Gene Expression and FoxO signaling pathway. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is PLK3.

抗原:Synthesized peptide derived from the Internal region of human Snk

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Expressed at higher level in the fetal lung, kidney, spleen and heart.