

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SPIC Polyclonal Antibody

多克隆抗体

产品货号: PA05896

第三版

描述: Controls the development of red pulp macrophages required for red blood cells recycling and iron homeostasis. Transcription factor that binds to the PU-box, a purine-rich DNA sequence (5'-GAGGA[AT]-3') that can act as a lymphoid-specific enhancer. Regulates VCAM1 gene expression (By similarity).

抗原: Synthesized peptide derived from the Internal region of human Spi-C

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Preferentially detected in fetal and adult spleen, lymph nodes and at lower levels in bone marrow and fetal liver.