

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SUPT3H Polyclonal Antibody

多克隆抗体

产品货号: PA05899

第三版

描述:SUPT3H (SPT3 Homolog, SAGA And STAGA Complex Component) is a Protein Coding gene. Diseases associated with SUPT3H include Blastic Plasmacytoid Dendritic Cell. Among its related pathways are Validated targets of C-MYC transcriptional activation and C-MYC pathway. GO annotations related to this gene include protein heterodimerization activity and histone acetyltransferase activity.

抗原:Synthesized peptide derived from the Internal region of human SPT3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in all tissues tested including pancreas, kidney, skeletal muscle, liver, lung, placenta, brain and heart.