

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey HTR1F Polyclonal Antibody

多克隆抗体

产品货号: PA05902

第三版

描述:HTR1F (5-Hydroxytryptamine Receptor 1F) is a Protein Coding gene. Diseases associated with HTR1F include Migraine With Or Without Aura 1. Among its related pathways are Serotonergic synapse and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and serotonin binding. An important paralog of this gene is HTR1E.

抗原:Synthesized peptide derived from the Internal region of human SR-1F

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: