

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STAT3 Polyclonal Antibody

多克隆抗体

产品货号: PA05928

第三版

描述: Transcription factor that binds to the interleukin-6 (IL-6)-responsive elements identified in the promoters of various acute-phase protein genes. Activated by IL31 through IL31RA.

抗原: Synthesized peptide derived from human Stat3 around the non-phosphorylation site of Tyr705.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.