

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STAT3 Polyclonal Antibody

多克隆抗体

产品货号: PA05929

第三版

**描述:** Transcription factor that binds to the interleukin-6 (IL-6)-responsive elements identified in the promoters of various acute-phase protein genes. Activated by IL31 through IL31RA.

**抗原:** Synthesized peptide derived from human Stat3 around the non-phosphorylation site of Ser727.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.