



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human STEAP1 Polyclonal Antibody

多克隆抗体

产品货号: PA05938

第三版

描述: Metalloreductase that has the ability to reduce both Fe(3+) to Fe(2+) and Cu(2+) to Cu(1+). Uses NAD(+) as acceptor.

抗原: Synthesized peptide derived from the N-terminal region of human STEAP1.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Ubiquitously expressed. Highly expressed in prostate tumors.