

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ELP2 Polyclonal Antibody

多克隆抗体

产品货号: PA05939

第三版

描述:Regulates the ligand-dependent activation of STAT3. Acts as subunit of the RNA polymerase II elongator complex, which is a histone acetyltransferase component of the RNA polymerase II (Pol II) holoenzyme and is involved in transcriptional elongation. Elongator may play a role in chromatin remodeling and is involved in acetylation of histones H3 and probably H4.

抗原:Synthesized peptide derived from the C-terminal region of human Stlp1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: