

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZNF280A Polyclonal Antibody

多克隆抗体

产品货号: PA05944

第三版

描述:This gene encodes a zinc finger protein. The encoded protein contains 4 C2H2-type zinc fingers, which are commonly found in transcription factors. A variety of functions may be performed by this type of zinc finger protein, including the binding of DNA or RNA. ZNF280A (Zinc Finger Protein 280A) is a Protein Coding gene. An important paralog of this gene is ZNF280B.

抗原:Synthesized peptide derived from the C-terminal region of human SUHW1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: