

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SFSWAP Polyclonal Antibody

多克隆抗体

产品货号: PA05947

第三版

描述: This gene encodes a human homolog of Drosophila splicing regulatory protein. This gene autoregulates its expression by control of splicing of its first two introns. In addition, it also regulates the splicing of fibronectin and CD45 genes. Two transcript variants encoding different isoforms have been identified. SFSWAP (Splicing Factor SWAP Homolog) is a Protein Coding gene. Among its related pathways are mRNA Splicing - Major Pathway. GO annotations related to this gene include RNA binding.

抗原: Synthesized peptide derived from the Internal region of human SWAP

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: