

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SS18 Polyclonal Antibody

多克隆抗体

产品货号: PA05961

第三版

描述:SS18 (SS18, NBAF Chromatin Remodeling Complex Subunit) is a Protein Coding gene. Diseases associated with SS18 include Sarcoma, Synovial and Sarcoma. Among its related pathways are Transcriptional misregulation in cancer and Chromatin Regulation / Acetylation. GO annotations related to this gene include ligand-dependent nuclear receptor transcription coactivator activity. An important paralog of this gene is SS18L1.

抗原:Synthesized peptide derived from the N-terminal region of human SYT

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Fairly ubiquitously expressed. Expressed in synovial sarcomas and in other human cell lines. The fusion genes SSXT-SSX1 and SSXT-SSX2 are expressed only in synovial sarcomas.