

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TAAR3P Polyclonal Antibody

多克隆抗体

产品货号: PA05965

第三版

描述:TAAR3P (Trace Amine Associated Receptor 3, Pseudogene) is a Pseudogene. Among its related pathways are Peptide ligand-binding receptors and Signaling by GPCR.

抗原:Synthesized peptide derived from the Internal region of human TAAR3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Not expressed in the pons, thalamus, globus pallidus, caudate, putamen or cerebellum.