

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TNKS1BP1 Polyclonal Antibody

多克隆抗体

产品货号: PA05966

第三版

**描述:** TNKS1BP1 (Tankyrase 1 Binding Protein 1) is a Protein Coding gene. Among its related pathways are Gene Expression and Deadenylation-dependent mRNA decay. GO annotations related to this gene include enzyme binding and ankyrin binding.

**抗原:** Synthesized peptide derived from the C-terminal region of human TAB182

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Detected in testis, ovary, lung, skeletal muscle, heart, prostate and pancreas, and at very low levels in brain and peripheral blood leukocytes.