

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TBC1D3A/B/C Polyclonal Antibody

多克隆抗体

产品货号: PA05981

第三版

描述: Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.

抗原: Synthesized peptide derived from the C-terminal region of human TBC1D3A/B/C

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in liver, skeletal muscle, kidney, pancreas, spleen, testis, ovary, small intestine and peripheral blood leukocytes. Overexpressed in prostate cancers.