

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TBPL2 Polyclonal Antibody

多克隆抗体

产品货号: PA05982

第三版

描述: TBPL2 (TATA-Box Binding Protein Like 2) is a Protein Coding gene. Among its related pathways are Respiratory electron transport, ATP synthesis by chemiosmotic coupling, and heat production by uncoupling proteins. and Influenza A. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is TBP.

抗原: Synthesized peptide derived from the Internal region of human TBPL2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitously expressed in all tissues examined with highest levels in heart, lung, ovary, spleen and testes.