

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TBX10 Polyclonal Antibody

多克隆抗体

产品货号: PA05984

第三版

**描述:** This gene encodes a member of the T-box family of transcription factors. These transcription factors share a DNA-binding domain called the T-box, and play a role in several developmental processes including early embryonic cell fate and organogenesis. The encoded protein is a member of the T-box 1 subfamily. Mutations in this gene are thought to be a cause of isolated cleft lip with or without cleft palate.

**抗原:** Synthesized peptide derived from the Internal region of human TBX10

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**