

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TBX15/18 Polyclonal Antibody

多克隆抗体

产品货号: PA05985

第三版

描述: TBX15 (T-Box 15) is a Protein Coding gene. Diseases associated with TBX15 include Cousin Syndrome and Acromegaloid Facial Appearance Syndrome. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is TBX18. TBX18 (T-Box 18) is a Protein Coding gene. Diseases associated with TBX18 include Pelviureteric Junction Obstruction and Huntington Disease. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is TBX15.

抗原: Synthesized peptide derived from the Internal region of human TBX15/18

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: