

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PHLDA1 Polyclonal Antibody

多克隆抗体

产品货号: PA05996

第三版

描述: This gene encodes an evolutionarily conserved proline-histidine rich nuclear protein. The encoded protein may play an important role in the anti-apoptotic effects of insulin-like growth factor-1. PHLDA1 (Pleckstrin Homology Like Domain Family A Member 1) is a Protein Coding gene. Diseases associated with PHLDA1 include Morpheiform Basal Cell Carcinoma and Juvenile Astrocytoma. Among its related pathways are Regulation of PLK1 Activity at G2/M Transition and Cell Cycle, Mitotic. An important paralog of this gene is PHLDA3.

抗原: Synthesized peptide derived from the C-terminal region of human TDAG51

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed with highest levels in pancreas. Strongly expressed by benign melanocytic nevi, and progressively reduced expressed in primary and metastatic melanomas (at protein level).