

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TNS1 Polyclonal Antibody

多克隆抗体

产品货号: PA06004

第三版

描述:The protein encoded by this gene localizes to focal adhesions, regions of the plasma membrane where the cell attaches to the extracellular matrix. This protein crosslinks actin filaments and contains a Src homology 2 (SH2) domain, which is often found in molecules involved in signal transduction. This protein is a substrate of calpain II. Alternative splicing results in multiple transcript variants encoding different isoforms. TNS1 (Tensin 1) is a Protein Coding gene. Diseases associated with TNS1 include Bannayan-Riley-Ruvalcaba Syndrome and Ruvalcaba Syndrome. Among its related pathways are Focal Adhesion and Integrin Pathway. GO annotations related to this gene include poly(A) RNA binding and actin binding. An important paralog of this gene is TNS3.

抗原:Synthesized peptide derived from human Tensin-1 around the non-phosphorylation site of Tyr1326.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.