

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TFEB Polyclonal Antibody

多克隆抗体

产品货号: PA06011

第三版

描述:TFEB (Transcription Factor EB) is a Protein Coding gene. Diseases associated with TFEB include Renal Cell Carcinoma, Papillary and Microphthalmia. Among its related pathways are Apoptosis and Autophagy and PI3K-AKT-mTOR signaling pathway and therapeutic opportunities. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein dimerization activity. An important paralog of this gene is MITF.

抗原:Synthesized peptide derived from the N-terminal region of human TFEB

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: