

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GTF3C3 Polyclonal Antibody

多克隆抗体

产品货号: PA06015

第三版

描述:The protein encoded by this gene is part of the TFIIC2 complex, which binds to the promoters of small nuclear and cytoplasmic RNA genes in order to recruit RNA polymerase III. The TFIIC2 complex is composed of six subunits. Two transcript variants encoding different isoforms have been found for this gene. GTF3C3 (General Transcription Factor IIIC Subunit 3) is a Protein Coding gene. Diseases associated with GTF3C3 include Tooth Agenesis. Among its related pathways are Gene Expression and RNA Polymerase III Transcription Initiation.

抗原:Synthesized peptide derived from the N-terminal region of human TFIIC102

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: