

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OXER1 Polyclonal Antibody

多克隆抗体

产品货号: PA06016

第三版

描述:OXER1 (Oxoeicosanoid Receptor 1) is a Protein Coding gene. Among its related pathways are Peptide ligand-binding receptors and Signaling by GPCR. GO annotations related to this gene include G-protein coupled receptor activity and 5-hydroxy-6E,8Z,11Z,14Z-icosatetraenoic acid binding. An important paralog of this gene is HCAR1.

抗原:Synthesized peptide derived from the Internal region of human TG1019

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in various tissues except brain. Expression is more intense in liver, kidney, peripheral leukocyte, lung, and spleen than in other tissues. Highly expressed in eosinophils, neutrophils, and lung macrophages.