

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MASTL Polyclonal Antibody

多克隆抗体

产品货号: PA06022

第三版

描述: This gene encodes a microtubule-associated serine/threonine kinase. Mutations at this locus have been associated with autosomal dominant thrombocytopenia, also known as thrombocytopenia-2. Alternatively spliced transcript variants have been described for this locus. MASTL (Microtubule Associated Serine/Threonine Kinase Like) is a Protein Coding gene. Diseases associated with MASTL include Thrombocytopenia 2 and Autosomal Thrombocytopenia With Normal Platelets. Among its related pathways are Mitotic Prophase and Cell Cycle, Mitotic. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is MAST4.

抗原: Synthesized peptide derived from the C-terminal region of human THC2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: