

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TIGD3 Polyclonal Antibody

多克隆抗体

产品货号: PA06027

第三版

描述:The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known. TIGD3 (Tigger Transposable Element Derived 3) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is TIGD4.

抗原:Synthesized peptide derived from the C-terminal region of human TIGD3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: